




STEAM InCator®


Class 4 Steam Chemical Indicator 9001



-  WAP - advanced waterproof technology, stable color, easy to read
-  Point-break strip design, all-round monitoring without blind spot
-  Eco-friendly lead-free formula to eliminate lead contamination


ISO 13485

Safer Supply, More Than Cleaner




InCator®

Class 4 Steam Chemical Integrator



Recommended Placement



Wrapped Towel Pack

Place a class 5 integrator in geometric center, not on top.



Rigid Container

Place two class 5 integrators in opposite corners of each level.



Peel Pouches

When an internal class 5 integrator is visible, an external type 1 CI is not required.



Non-woven Wrap

Place a class 5 integrator in geometric center, not on top.

Ordering information

Model	Product Name	Size(inchs)	PCS(Bag)	Bags(Case)
9001	InCator® Class 4 Steam Chemical Indicator	7.87*0.63	500pcs/bag	10bags/case

InCator®

Class 4 Steam Chemical Integrator

JUSHA Steam Chemical Indicator Strip 9001

Unprocessed



PASS



FALL



InCator[®] Class 4 Steam Chemical Integrator

JUSHA Steam Chemical Indicator Strip 9001

Unprocessed



JUSHA **Steam Chemical Indicator**
InCator[®] 9001 Strip turns dark when processing

STEAM

JUSHA **Steam Chemical Indicator**
InCator[®] 9001 Strip turns dark when processing

PASS



JUSHA **Steam Chemical Indicator**
InCator[®] 9001 Strip turns dark when processing

STEAM

JUSHA **Steam Chemical Indicator**
InCator[®] 9001 Strip turns dark when processing

FALL



JUSHA **Steam Chemical Indicator**
InCator[®] 9001 Strip turns dark when processing

STEAM

JUSHA **Steam Chemical Indicator**
InCator[®] 9001 Strip turns dark when processing